

For Research Use Only

Envoplakin Polyclonal antibody

Catalog Number: 25441-1-AP



Basic Information

Catalog Number:

25441-1-AP

Size:

150ul, Concentration: 350 ug/ml by Nanodrop and 287 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22188

GenBank Accession Number:

BC126103

GeneID (NCBI):

2125

UNIPROT ID:

Q92817

Full Name:

envoplakin

Calculated MW:

2033 aa, 232 kDa

Observed MW:

210 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A431 cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 25441-1-AP (Envoplakin antibody at dilution of 1:300 incubated at room temperature for 1.5 hours.