

For Research Use Only

SSTR5 Polyclonal antibody

Catalog Number: 51029-2-AP



Basic Information

Catalog Number: 51029-2-AP	GenBank Accession Number: BC146576	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 347 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 6755	Recommended Dilutions: WB 1:500-1:2400
Source: Rabbit	Full Name: somatostatin receptor 5	
Isotype: IgG	Calculated MW: 364 aa, 39 kDa	
Immunogen Catalog Number: AG0519	Observed MW: 50 kDa, 65 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human placenta tissue,
Species Specificity: human, mouse	

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 51029-2-AP (SSTR5 antibody) at dilution of 1:0 incubated at room temperature for 1.5 hours.