

For Research Use Only

SSTR5 Polyclonal antibody

Catalog Number: 51029-2-AP



Basic Information

Catalog Number:

51029-2-AP

Size:

150ul , Concentration: 347 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0519

GenBank Accession Number:

BC146576

GeneID (NCBI):

6755

UNIPROT ID:

P35346

Full Name:

somatostatin receptor 5

Calculated MW:

364 aa, 39 kDa

Observed MW:

50 kDa, 65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2400

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : human placenta tissue,

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 51029-2-AP (SSTR5 antibody) at dilution of 1:0 incubated at room temperature for 1.5 hours.