

For Research Use Only

Jph2 Polyclonal antibody

Catalog Number: 51054-2-AP



Basic Information

Catalog Number: 51054-2-AP	GenBank Accession Number: NM_021566	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 433 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 59091	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	Full Name: junctophilin 2	
Isotype: IgG	Calculated MW: 75 kDa 45 kDa	
Immunogen Catalog Number: AG0468		

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human pancreas cancer tissue,
Species Specificity: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human pancreas cancer using 51054-2-AP (Jph2 antibody) at dilution of 1:50 (under 10x lens).