

For Research Use Only

# PE Anti-Human CD1c (L161)

Catalog Number: PE-65176



## Basic Information

**Catalog Number:**  
PE-65176

**Size:**  
100tests, 10 µl/test

**Source:**  
Mouse

**Isotype:**  
Mouse IgG1 kappa

**GenBank Accession Number:**  
BC126465

**GeneID (NCBI):**  
911

**Full Name:**  
CD1c molecule

**Purification Method:**

The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

**CloneNo.:**  
L161

**Excitation/Emission maxima wavelengths:**  
496 nm, 565 nm / 578 nm

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Human

## Background Information

### Storage


**Storage:**  
Store at 2-8°C. Avoid exposure to light.

**Storage Buffer:**  
PBS with 0.1% sodium azide.


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone:HiB19) and 10.00 ul/test PE-Anti-Human CD1c (PE-65176, Clone:L161). Cells were not fixed.



1X10<sup>6</sup> human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone:HiB19) and PE-mouse IgG1 isotype control antibody. Cells were not fixed.