



FlexAble^{2.0}

Antibody Labeling Kits

Any antibody. Any color. Any time.

Quick and easy antibody labeling with even more dyes
for brighter signal in fluorescence applications.

Background image: Mouse colon: anti-FUS (11570-1-AP) + FlexAble 2.0 CoraLite® Plus 647 Kit (KFA503, magenta), anti- α -SMA (67735-1-Ig) + FlexAble 2.0 CoraLite® Plus 555 Kit (KFA522, orange), anti-E-cadherin (60335-1-Ig) + FlexAble 2.0 CoraLite® Plus 488 Kit (KFA561, green), and DAPI (blue).

ptglab.com

FlexAble^{2.0}

Antibody Labeling Kits

Any antibody. Any color. Any time.

- No need for antibody purification before or after labeling
- Site-specific antibody labeling with consistent DOL
- Label as little as 0.5ug and scale up to fit your needs
- Perfect tool for multiplexing, even within the same IgG species

COMPATIBLE WITH:

- ✓ Antibodies from any supplier
- ✓ Any antibody concentration
- ✓ Antibodies that contain BSA and glycerol

SPECIFICITIES

Rabbit IgG
Mouse IgG1
Mouse IgG2a
Mouse IgG2b
Human
Rat Kappa LC

CONJUGATES

CoraLite[®] Plus 405
CoraLite[®] Plus 488
CoraLite[®] Plus 555
CoraLite[®] Plus 594
CoraLite[®] Plus 647
CoraLite[®] Plus 750
FITC Plus

HRP
Biotin

