

Speedy ELISA Kit Offering

| Target | Catalog | Sensitivity |
|--------------------|---------|-------------|
| Human APOE | SE50006 | 0.03 ng/mL |
| Human HMGB1 | SE50009 | 0.05 ng/mL |
| Human AFP | SE50010 | 7.32 ng/mL |
| Human TNF-alpha | SE50011 | 2.81 ng/mL |
| Human OX40L/TNFSF4 | SE50012 | 0.03 ng/mL |
| Human IL-6 | SE50013 | 1.06 ng/mL |
| Mouse IL-2 | SE51002 | 1.71 ng/mL |
| Rat IL-2 | SE52002 | 2.67 ng/mL |

Speedy Kits enable analyte detection in 90 minutes or less, nearly three times faster than conventional ELISAs. Simply combine samples, standards, and antibodies in a single incubation step, followed by washing, development, and reading.

