



FlexAble^{2.0}

Antibody Labeling Kits

Any antibody. Any color. Any time.

Quick and easy antibody labeling with even more dyes for brighter signals for immunofluorescence, western blot, and flow cytometry. Conjugate fluorochromes, enzymes, and molecules in any buffer condition.

- ✓ Fast and easy 2-step protocol
- ✓ Ready to use in 10 minutes
- ✓ No buffer exchange required

SPECIFICITIES

Rabbit IgG
Mouse IgG1
Mouse IgG2a
Mouse IgG2b
Human
Rat Kappa LC

CONJUGATES

CoraLite[®] Plus 405
CoraLite[®] Plus 488
CoraLite[®] Plus 555
CoraLite[®] Plus 594
CoraLite[®] Plus 647
CoraLite[®] Plus 750
FITC Plus

HRP
Biotin

