

For Research Use Only

ICP27 Polyclonal antibody

Catalog Number: 10186-1-AP **1 Publications**



Basic Information

Catalog Number: 10186-1-AP
Size: 150ul , Concentration: 800 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard;
Source: Rabbit
Isotype: IgG
Immunogen Catalog Number: AG0239

GenBank Accession Number: AAF43147
GeneID (NCBI):
Full Name:
Calculated MW: 55 kDa

Purification Method: Antigen affinity purification

Applications

Tested Applications: ELISA
Cited Applications: WB
Species Specificity: hsv

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
M L Li	23237087	Acta Virol	WB

Storage

Storage: Store at -20°C. Stable for one year after shipment.
Storage Buffer: PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data