

For Research Use Only

Goat-IgG Polyclonal antibody

Catalog Number: 10285-1-AP

48 Publications



Basic Information

Catalog Number:

10285-1-AP

GenBank Accession Number:

GeneID (NCBI):

Purification Method:

Antigen affinity purification

Size:

150ul , Concentration: 12000 µg/ml
by Nanodrop and 400 µg/ml by
Bradford method using BSA as the
standard;

Full Name:

Calculated MW:

155 kDa

Source:

Rabbit

Isotype:

IgG

Applications

Tested Applications:

ELISA

Cited Applications:

IF, RIP, WB

Species Specificity:

goat

Cited Species:

human, rat, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jian Wen	32996777	Cancer Biother Radiopharm	WB
Lan Zhang	34586738	Clin Transl Med	RIP
Hucheng Chen	30272275	Mol Med Rep	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

