

For Research Use Only

# BL2 1-lysate Polyclonal antibody

Catalog Number: 10557-1-Ig



## Basic Information

<b>Catalog Number:</b> 10557-1-Ig	<b>GenBank Accession Number:</b> GeneID (NCBI):	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 150ul , Concentration: 1800 µg/ml by Nanodrop and 700 µg/ml by Bradford method using BSA as the standard;	<b>Full Name:</b> <b>Calculated MW:</b> 1-300 kDa	<b>Recommended Dilutions:</b> WB 1:1000-1:8000
<b>Source:</b> Rabbit		
<b>Isotype:</b> IgG		
<b>Immunogen Catalog Number:</b> AG0049		

## Applications

<b>Tested Applications:</b> WB,ELISA	<b>Positive Controls:</b> WB : Total bacteria lysates of BL21,
<b>Species Specificity:</b> bacteria	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data