

For Research Use Only

ARL6IP4 Polyclonal antibody

Catalog Number: 10563-1-AP



Basic Information

Catalog Number: 10563-1-AP	GenBank Accession Number: BC001958	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 587 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 51329	
Source: Rabbit	Full Name: ADP-ribosylation-like factor 6 interacting protein 4	
Isotype: IgG	Calculated MW: 25 kDa	
Immunogen Catalog Number: AG0863		

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

ARL6IP4 is recognized to be a binding protein of ARL6 that is another member in the Ras family, implicating a possibility of mediating Ras signaling. ARL6 gene mutation has been linked to Bardet-Biedl syndrome type 3 which is a multisystemic disorder characterized by defects in the microtubule-based transport process.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data