

For Research Use Only

# ERAL1 Polyclonal antibody

Catalog Number: 11478-1-AP

Featured Product

12 Publications



## Basic Information

### Catalog Number:

11478-1-AP

### Size:

150ul, Concentration: 347 µg/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG2019

### GenBank Accession Number:

BC019094

### GeneID (NCBI):

26284

### Full Name:

Era G-protein-like 1 (E. coli)

### Calculated MW:

437 aa, 48 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

ELISA

### Cited Applications:

IF, WB

### Species Specificity:

human

### Cited Species:

human, mouse

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Agata Rozanska	28512204	Biochem J	WB
Tian Luo	28553621	Front Cell Infect Microbiol	WB, IF
Di Hu	30877431	Acta Neuropathol	

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

