

For Research Use Only

# FGF19 Polyclonal antibody

Catalog Number: 12291-1-AP



## Basic Information

Catalog Number:  
12291-1-AP

Size:  
150ul , Concentration: 260 µg/ml by  
Bradford method using BSA as the  
standard;

Source:  
Rabbit

Isotype:  
IgG

Immunogen Catalog Number:  
AG2978

GenBank Accession Number:  
BC017664

GeneID (NCBI):  
9965

Full Name:  
fibroblast growth factor 19

Calculated MW:  
216 aa, 24 kDa

Purification Method:  
Antigen affinity purification

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data