

For Research Use Only

# FGF13 Polyclonal antibody

Catalog Number: 13201-1-AP

3 Publications



## Basic Information

<b>Catalog Number:</b> 13201-1-AP	<b>GenBank Accession Number:</b> BC034340	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 2258	
<b>Source:</b> Rabbit	<b>Full Name:</b> fibroblast growth factor 13	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 28 kDa	
<b>Immunogen Catalog Number:</b> AG3942		

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
IF, WB

**Species Specificity:**  
human, mouse, rat

**Cited Species:**  
human, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhang Xiujun X	23554914	PLoS One	WB
Xiaojing Yue	29405966	Placenta	WB,IF
J-J Song	30720165	Eur Rev Med Pharmacol Sci	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Aliquoting is unnecessary for -20°C storage**

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

