

For Research Use Only

# GDF5 Polyclonal antibody

Catalog Number: 14645-1-AP **1 Publications**



## Basic Information

|  |  |  |
|--|--|--|
| <b>Catalog Number:</b><br>14645-1-AP   | <b>GenBank Accession Number:</b><br>BC032495         | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>150ul , Concentration: 200 µg/ml by Bradford method using BSA as the standard; | <b>GeneID (NCBI):</b><br>8200                        |  |
| <b>Source:</b><br>Rabbit   | <b>Full Name:</b><br>growth differentiation factor 5 |  |
| <b>Isotype:</b><br>IgG   | <b>Calculated MW:</b><br>55 kDa                      |  |
| <b>Immunogen Catalog Number:</b><br>AG6238   |  |  |

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
IHC

**Species Specificity:**  
human, mouse, rat

**Cited Species:**  
human

## Background Information

### Notable Publications

| Author                | Pubmed ID | Journal                  | Application |
|-----------------------|-----------|--------------------------|-------------|
| Phaedra R Theodoridis | 33568529  | Proc Natl Acad Sci U S A | IHC         |

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data