

For Research Use Only

# DYNLRB1 Polyclonal antibody

Catalog Number: 15049-1-AP

3 Publications



## Basic Information

Catalog Number:

15049-1-AP

Size:

150ul , Concentration: 200 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7032

GenBank Accession Number:

BC002481

GeneID (NCBI):

83658

Full Name:

dynein, light chain, roadblock-type 1

Calculated MW:

11 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

IF

Species Specificity:

human, mouse, rat

Cited Species:

human, rat

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Syeda Ridita Sharif	27646688	Mol Cells	IF
Md Ariful Islam	26467288	Mol Cells	IF
Md Ariful Islam	26272270	Exp Mol Med	IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

