

For Research Use Only

# FGF6-Specific Polyclonal antibody

Catalog Number: 16878-1-AP



## Basic Information

|                               |  |   |
|-------------------------------|--|---|
| Catalog Number:<br>16878-1-AP | GenBank Accession Number:<br>NM_020996   | Purification Method:<br>Antigen affinity purification |
| Size:<br>150ul                | GeneID (NCBI):<br>2251                   |   |
| Source:<br>Rabbit             | Full Name:<br>fibroblast growth factor 6 |   |
| Isotype:<br>IgG               | Calculated MW:<br>23 kDa                 |   |

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human, mouse, rat

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data