

For Research Use Only

GAL3ST1 Polyclonal antibody

Catalog Number: 18503-1-AP



Basic Information

Catalog Number: 18503-1-AP	GenBank Accession Number: BC014649	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 240 µg/ml by Nanodrop and 200 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 9514	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: galactose-3-O-sulfotransferase 1	
Isotype: IgG	Calculated MW: 49 kDa	
Immunogen Catalog Number: AG13353	Observed MW: 48 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, HEK-293 cells, SGC-7901 cells
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 18503-1-AP (GAL3ST1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.