

For Research Use Only

ENTPD4 Polyclonal antibody

Catalog Number: 19597-1-AP



Basic Information

Catalog Number: 19597-1-AP	GenBank Accession Number: BC034477	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Nanodrop and 287 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 9583	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: ectonucleoside triphosphate diphosphohydrolase 4	
Isotype: IgG	Calculated MW: 70 kDa	
Immunogen Catalog Number: AG5175	Observed MW: 70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells, human placenta tissue
Species Specificity: human, mouse	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 19597-1-AP (ENTPD4 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.