

For Research Use Only

IGF2 Polyclonal antibody

Catalog Number: 24911-1-AP **1 Publications**



Basic Information

Catalog Number: 24911-1-AP	GenBank Accession Number: BC000531	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 280 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 3481	
Source: Rabbit	Full Name: insulin-like growth factor 2 (somatomedin A)	
Isotype: IgG	Calculated MW: 236 aa, 26 kDa	
Immunogen Catalog Number: AG21397		

Applications

Tested Applications:
IHC, ELISA

Cited Applications:
IHC

Species Specificity:
human

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jinglin Zhang	36373776	J Pathol	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data