

For Research Use Only

# SPATC1 Polyclonal antibody

Catalog Number: 25861-1-AP **2 Publications**



## Basic Information

<b>Catalog Number:</b> 25861-1-AP	<b>GenBank Accession Number:</b> BC053547	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 550 µg/ml by Nanodrop and 307 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 375686	<b>Recommended Dilutions:</b> WB 1:500-1:1000
<b>Source:</b> Rabbit	<b>Full Name:</b> spermatogenesis and centriole associated 1	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 500 aa, 52 kDa	
<b>Immunogen Catalog Number:</b> AG12157	<b>Observed MW:</b> 60 kDa	

## Applications

<b>Tested Applications:</b> WB,ELISA	<b>Positive Controls:</b> WB : human testis tissue,
<b>Cited Applications:</b> IF	
<b>Species Specificity:</b> human, mouse	
<b>Cited Species:</b> mouse	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Brett Nixon	33451771	Mol Cell Proteomics	IF
Brett Nixon	30072580	Mol Cell Proteomics	IF

## Storage


**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 25861-1-AP (SPATC1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.