

For Research Use Only

TTL10 Polyclonal antibody

Catalog Number: 26014-1-AP



Basic Information

Catalog Number: 26014-1-AP	GenBank Accession Number: BC126152	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 µg/ml by Nanodrop and 233 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 254173	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: tubulin tyrosine ligase-like family, member 10	
Isotype: IgG	Observed MW: 45 kDa	
Immunogen Catalog Number: AG23258		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, rat testis tissue
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 26014-1-AP (TLL10 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.