

For Research Use Only

ACOT7L Polyclonal antibody

Catalog Number: 51016-2-AP



Basic Information

Catalog Number: 51016-2-AP	GenBank Accession Number: AK125299	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop;	GeneID (NCBI): 344967	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: acyl-CoA thioesterase 7 pseudogene	
Isotype: IgG	Calculated MW: 28 kDa	
	Observed MW: 28 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, human brain tissue, Hela cells
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 51016-2-AP (ACOT7L antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.