

For Research Use Only

Prnd Polyclonal antibody

Catalog Number: 51026-2-AP



Basic Information

Catalog Number: 51026-2-AP	GenBank Accession Number: BC025140	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 26434	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: prion protein dublet	
Isotype: IgG	Calculated MW: 21 kDa	
Immunogen Catalog Number: AG0731	Observed MW: 34 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : human skeletal muscle tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 51026-2-AP (Prnd antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.