

For Research Use Only

TWEAKR Monoclonal antibody

Catalog Number: 60104-1-Ig



Basic Information

Catalog Number: 60104-1-Ig	GenBank Accession Number: BC002718	Purification Method: Caprylic Acid/Ammonium Sulfate Precipitation
Size: 150ul	GeneID (NCBI): 51330	CloneNo.: 2E10C10
Source: Mouse	Full Name: tumor necrosis factor receptor superfamily, member 12A	
Isotype: IgG2b	Calculated MW: 129 aa, 14 kDa	
Immunogen Catalog Number: AG2689		

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human heart tissue were subjected to SDS PAGE followed by western blot with 60104-1-Ig (TWEAKR antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.