

For Research Use Only

DR6 Monoclonal antibody

Catalog Number: 60107-1-Ig



Basic Information

Catalog Number: 60107-1-Ig	GenBank Accession Number: BC005192	Purification Method: Caprylic acid/ammonium sulfate precipitation
Size: 150ul , Concentration: 600 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 27242	CloneNo.: 4D4F8
Source: Mouse	Full Name: tumor necrosis factor receptor superfamily, member 21	
Isotype: IgG1	Calculated MW: 15 kDa	
Immunogen Catalog Number: AG1126	Observed MW: 72 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data