



IHCeasy SAMD4B Ready-To-Use IHC Kit

Catalog Number: KHC1776

General Information

Sample type: FFPE tissue Cited sample type: Reactivity: Human, Mouse, Rat Cited Reactivity: Assay type: Immunohistochemistry Primary antibody type: Rabbit Polyclonal

Secondary antibody type: Polymer-HRP-Goat anti-Rabbit

Kit Component

Antigen Retrieval Buffer 100 mL	50×
Washing Buffer 100 mL × 2	20×
Blocking Buffer 5 mL	RTU
Primary Antibody 5 mL	RTU
Secondary Antibody 5 mL	RTU
Chromogen Component A 0.2 mL	RTU
Chromogen Component B 4 mL	RTU
Signal Enhancer 5 mL	RTU
Counter Staining Reagent 5 mL	RTU
Mounting Media 5 mL	RTU
Control Slide 1 slide (Optional)	FFPE
Datasheet 1 Copy	
Manual 1 Copy	


Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Synonyms

hSmaug2, Protein Smaug homolog 2, SAMD4B, Smaug 2, SMAUG2, SMGB


Selected Validation Data




Immunohistochemical analysis of paraffinembedded human thyroid cancer tissue slide using KHC1776 (SAMD4B IHC Kit).




Immunohistochemical analysis of paraffinembedded human urothelial carcinoma tissue slide using KHC1776 (SAMD4B IHC Kit).




Immunohistochemical analysis of paraffinembedded mouse small intestine tissue slide using KHC1776 (SAMD4B IHC Kit).



Immunohistochemical analysis of paraffinembedded mouse kidney tissue slide using KHC1776 (SAMD4B IHC Kit).



Immunohistochemical analysis of paraffinembedded mouse pancreas tissue slide using KHC1776 (SAMD4B IHC Kit).



Immunohistochemical analysis of paraffinembedded rat small intestine tissue slide using KHC1776 (SAMD4B IHC Kit).