

For Research Use Only

BCRP/ABCG2 Polyclonal antibody



Catalog Number: 10051-1-AP **8 Publications**

Basic Information

Catalog Number: 10051-1-AP
Size: 150UL
Source: Rabbit
Isotype: IgG
Immunogen Catalog Number: AG0005

GenBank Accession Number: BC021281
GeneID (NCBI): 9429
Full Name: ATP-binding cassette, sub-family G (WHITE), member 2
Calculated MW: 70 kDa

Purification Method: Antigen affinity purification

Applications

Tested Applications: ELISA
Cited Applications: WB
Species Specificity: human
Cited Species: canine, human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jingsheng Yuan	31078735	Cancer Lett	WB
Manxin Xie	29567542	Biomed Pharmacother	WB
Ping Xie	28569775	Cell Death Dis	WB

Storage

Storage: Store at -20°C. Stable for one year after shipment.
Storage Buffer: PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data