

For Research Use Only

RTN4/NOGO Polyclonal ANTIBODY

Catalog Number: 10367-1-AP



Basic Information

Catalog Number:
10367-1-AP

Size:
150UL, Concentration: 133 µg/ml by
Bradford method using BSA as the
standard;

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG0495

GenBank Accession Number:
BC001035

GeneID (NCBI):
57142

Full Name:
reticulon 4

Calculated MW:
22 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data