

For Research Use Only

# BBS7 Polyclonal ANTIBODY



Catalog Number: 13458-1-AP

## Basic Information

Catalog Number:

13458-1-AP

Size:

150UL, Concentration: 133 µg/ml by  
Bradford method using BSA as the  
standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4314

GenBank Accession Number:

BC032691

GeneID (NCBI):

55212

Full Name:

Bardet-Biedl syndrome 7

Calculated MW:

715 aa, 80 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF 1:20-1:200

## Applications

Tested Applications:

IF, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

IF : hTERT-RPE1 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



IF result (pericentriolar staining) of anti-BBS7 (13458-1-AP, 1:50) with serum-starved hTERT-RPE1 by Dr. Moshe Kim.