

For Research Use Only

# DNMT3A Polyclonal antibody

Catalog Number: 14623-1-AP

1 Publications



## Basic Information

Catalog Number:

14623-1-AP

Size:

150ul , Concentration: 127 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6271

GenBank Accession Number:

BC043617

GeneID (NCBI):

1788

UNIPROT ID:

Q9Y6K1

Full Name:

DNA (cytosine-5-)-methyltransferase

3 alpha

Calculated MW:

912 aa, 102 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

## Notable Publications

| Author   | Pubmed ID | Journal     | Application |
|----------|-----------|-------------|-------------|
| Yanna Ma | 35666085  | J Chemother | WB          |

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

