

For Research Use Only

# PRKAB1 Polyclonal ANTIBODY

Catalog Number: 15419-1-AP



## Basic Information

Catalog Number:  
15419-1-AP

Size:  
10 µg/150 µl

Source:  
Rabbit

Isotype:  
IgG

Purification Method:  
Antigen affinity purification

Immunogen Catalog Number:  
AG7707

GenBank Accession Number:  
BC001056

GeneID (NCBI):  
5564

Full Name:  
protein kinase, AMP-activated, beta 1 non-catalytic subunit

Calculated MW:  
30 kDa

Observed MW:

Recommended Dilutions:

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human,mouse,rat

**Note: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

## Background Information

## Notable Publications

| Author | Pubmed ID | Journal | Application |
|--------|-----------|---------|-------------|
|--------|-----------|---------|-------------|

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data