

For Research Use Only

ELF4 Polyclonal ANTIBODY

Catalog Number: 16020-1-AP



Basic Information

Catalog Number:
16020-1-AP

Size:
48 µg/150 µl

Source:
Rabbit

Isotype:
IgG

Purification Method:
Antigen affinity purification
Immunogen Catalog Number:
AG8836

GenBank Accession Number:
BC017194

GeneID (NCBI):
2000

Full Name:
E74-like factor 4 (ets domain transcription factor)

Calculated MW:
663aa, 71 kDa

Observed MW:

Recommended Dilutions:

Applications

Tested Applications:
ELISA

Species Specificity:
human,mouse,rat

Positive Controls:

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
--------	-----------	---------	-------------

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data