

For Research Use Only

ASB1 Polyclonal antibody

Catalog Number: 16095-1-AP

1 Publications



Basic Information

Catalog Number:

16095-1-AP

Size:

150ul, Concentration: 600 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8848

GenBank Accession Number:

BC014528

GeneID (NCBI):

51665

UNIPROT ID:

Q9Y576

Full Name:

ankyrin repeat and SOCS box-containing 1

Calculated MW:

335 aa, 37 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB: mouse spleen tissue, mouse testis tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Jinxing Lv	39733518	Redox Biol	WB, IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 16095-1-AP (ASB1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.