

For Research Use Only

HOXD4 Polyclonal antibody

Catalog Number: 17780-1-AP



Basic Information

Catalog Number: 17780-1-AP	GenBank Accession Number: BC016763	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Nanodrop and 213 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 3233	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: homeobox D4	
Isotype: IgG	Calculated MW: 28 kDa	
Immunogen Catalog Number: AG11134	Observed MW: 28 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : A549 cells,
Species Specificity: human	

Background Information


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 17780-1-AP (HOXD4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.