

For Research Use Only

WNT3 Polyclonal antibody

Catalog Number: 17983-1-AP **3 Publications**



Basic Information

Catalog Number: 17983-1-AP
Size: 150UL, Concentration: 160 µg/ml by Bradford method using BSA as the standard;
Source: Rabbit
Isotype: IgG
Immunogen Catalog Number: AG12399

GenBank Accession Number: BC114219
GeneID (NCBI): 7473
Full Name: wingless-type MMTV integration site family, member 3
Calculated MW: 355 aa, 40 kDa

Purification Method: Antigen affinity purification

Applications

Tested Applications: ELISA
Cited Applications: WB
Species Specificity: human, mouse
Cited Species: human, mouse, swine

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Qing Zuo	27853262	Sci Rep	WB
Min Li	27247798	J Ophthalmol	WB
Xiuying Zhong	30527743	Cell Metab	WB

Storage

Storage: Store at -20°C. Stable for one year after shipment.
Storage Buffer: PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data