

For Research Use Only

WNT3 Polyclonal ANTIBODY



Catalog Number: 17983-1-AP

3 Publications

Basic Information

Catalog Number:
17983-1-AP

Size:
24 µg/150 µl

Source:
Rabbit

Isotype:
IgG

Purification Method:
Antigen affinity purification

Immunogen Catalog Number:
AG12399

GenBank Accession Number:
BC114219

GeneID (NCBI):
7473

Full Name:
wingless-type MMTV integration site family,
member 3

Calculated MW:
355aa, 40 kDa

Observed MW:

Recommended Dilutions:

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse

Cited Species:
human, mouse, pig

Note: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Qing Zuo	27853262	Sci Rep	WB
Min Li	27247798	J Ophthalmol	WB
Xuying Zhong	30527743	Cell Metab	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data