

For Research Use Only

ADH1A Polyclonal ANTIBODY

Catalog Number: 17998-1-AP



Basic Information

Catalog Number:

17998-1-AP

Size:

150UL

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12467

GenBank Accession Number:

BC074738

GeneID (NCBI):

124

Full Name:

alcohol dehydrogenase 1A (class I),
alpha polypeptide

Calculated MW:

375 aa, 40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data