

For Research Use Only

Ephrin A1 Polyclonal ANTIBODY

Catalog Number: 19075-1-AP



Basic Information

Catalog Number:
19075-1-AP

Size:
0 µg/150 µl

Source:
Rabbit

Isotype:
IgG

Purification Method:
Antigen affinity purification
Immunogen Catalog Number:
AG4631

GenBank Accession Number:
BC032698

GeneID (NCBI):
1942

Full Name:
ephrin-A1

Calculated MW:
205aa, 24 kDa

Observed MW:

Recommended Dilutions:

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Note: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|--------|-----------|---------|-------------|
|--------|-----------|---------|-------------|

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data