

For Research Use Only

MGST3 Polyclonal ANTIBODY

Catalog Number:20114-1-AP



Basic Information

Catalog Number:
20114-1-AP

Size:
0 µg/150 µl

Source:
Rabbit

Isotype:
IgG

Purification Method:
Antigen affinity purification

Immunogen Catalog Number:
AG13932

GenBank Accession Number:
BC000505

GeneID (NCBI):
4259

Full Name:
microsomal glutathione S-transferase 3

Calculated MW:
152aa, 17 kDa

Observed MW:

Recommended Dilutions:

Applications

Tested Applications:
ELISA

Species Specificity:
human,mouse,rat

Note: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
--------	-----------	---------	-------------

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data