

For Research Use Only

# NELL1 Polyclonal antibody

Catalog Number: 21783-1-AP

3 Publications



## Basic Information

Catalog Number:

21783-1-AP

Size:

150ul, Concentration: 167 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16542

GenBank Accession Number:

BC069674

GeneID (NCBI):

4745

UNIPROT ID:

Q92832

Full Name:

NEL-like 1 (chicken)

Calculated MW:

810 aa, 90 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

## Notable Publications

Author	Pubmed ID	Journal	Application
Mengmeng Liu	26456257	J Endod	IHC
Mengyue Li	34512380	Front Physiol	
Tang Ronghua R	23264108	J Mol Histol	IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

