

For Research Use Only

HOXB4 Polyclonal antibody

Catalog Number: 22256-1-AP

1 Publications



Basic Information

Catalog Number:

22256-1-AP

Size:

150ul, Concentration: 650 ug/ml by Nanodrop and 280 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17761

GenBank Accession Number:

BC049204

GeneID (NCBI):

3214

UNIPROT ID:

P17483

Full Name:

homeobox B4

Calculated MW:

251 aa, 28 kDa

Observed MW:

28-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB: Jurkat cells,

Notable Publications

Author	Pubmed ID	Journal	Application
Shengtou Ye	40099691	Mol Carcinog	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 22256-1-AP (HOXB4 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.