

For Research Use Only

GADD45B Polyclonal antibody

Catalog Number: 22645-1-AP



Basic Information

Catalog Number: 22645-1-AP	GenBank Accession Number: BC113466	Purification Method: Antigen Affinity purified
Size: 150ul , Concentration: 340 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 4616	
Source: Rabbit	Full Name: growth arrest and DNA-damage-inducible, beta	
Isotype: IgG	Calculated MW: 160 aa, 18 kDa	
Immunogen Catalog Number: AG18524		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

