

For Research Use Only

KRT71 Polyclonal antibody

Catalog Number: 26149-1-AP



Basic Information

Catalog Number:	26149-1-AP	GenBank Accession Number:	BC103917	Purification Method:	Antigen affinity purification
Size:	150ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI):	112802	Recommended Dilutions:	WB 1:500-1:2000
Source:	Rabbit	UNIPROT ID:	Q3SY84		
Isotype:	IgG	Full Name:	keratin 71		
Immunogen Catalog Number:	AG23536	Observed MW:	57 kDa		

Applications

Tested Applications:	WB, ELISA	Positive Controls:	WB : mouse skin tissue,
Species Specificity:	human, mouse		

Background Information

KRT71 belongs to a family of type II keratins that are specifically expressed in the inner root sheath of hair follicles. The type II keratins are clustered in a region of chromosome 12q13.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 26149-1-AP (KRT71 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.