

For Research Use Only

ZRANB1 Polyclonal antibody

Catalog Number: 27050-1-AP



Basic Information

Catalog Number:

27050-1-AP

Size:

150ul , Concentration: 600 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25789

GenBank Accession Number:

AL832925

GeneID (NCBI):

54764

UNIPROT ID:

Q9UGIO

Full Name:

zinc finger, RAN-binding domain containing 1

Calculated MW:

81 kDa

Observed MW:

81 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human, mouse

Positive Controls:

IP : HeLa cells, mouse cerebellum tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-ZRANB1 (IP:27050-1-AP, 4ug; Detection:27050-1-AP 1:500) with HeLa cells lysate 1600 ug.



IP result of anti-ZRANB1 (IP:27050-1-AP, 4ug; Detection:27050-1-AP 1:500) with mouse cerebellum tissue lysate 2000 ug.