

For Research Use Only

TIAM1 Polyclonal antibody

Catalog Number: 27694-1-AP



Basic Information

Catalog Number:

27694-1-AP

GenBank Accession Number:

NM_003253

Purification Method:

Antigen affinity purification

Size:

150ul , Concentration: 1200 µg/ml by
Nanodrop;

GeneID (NCBI):

7074

Recommended Dilutions:

WB 1:500-1:1000

Source:

Rabbit

Full Name:

T-cell lymphoma invasion and
metastasis 1

Isotype:

IgG

Calculated MW:

178 kDa

Immunogen Catalog Number:

AG26777

Observed MW:

178 kDa

Applications

Tested Applications:

WB,ELISA

Positive Controls:

WB : HCT 116 cells, HEK-293T cells

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HCT 116 cells were subjected to SDS PAGE followed by western blot with 27694-1-AP (TIAM1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 27694-1-AP (TIAM1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.