

For Research Use Only

ABI3BP Polyclonal ANTIBODY

Catalog Number: 27727-1-AP



Basic Information

Catalog Number:

27727-1-AP

Size:

150UL, Concentration: 367 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26841

GenBank Accession Number:

BC030221

GeneID (NCBI):

25890

Full Name:

Abi family, member 3 (NESH) binding protein

Calculated MW:

119 kDa

Observed MW:

54 kDa, 72 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, Mouse

Positive Controls:

WB: C2C12 cell, mouse thymus tissue, mouse testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27727-1-AP (ABI3BP antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.