

For Research Use Only

KAT2B/PCAF Polyclonal antibody

Catalog Number: 28770-1-AP

6 Publications



Basic Information

Catalog Number:

28770-1-AP

Size:

150ul, Concentration: 700 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30238

GenBank Accession Number:

BC060823

GeneID (NCBI):

8850

UNIPROT ID:

Q92831

Full Name:

K(lysine) acetyltransferase 2B

Calculated MW:

93 kDa

Observed MW:

95-110 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB : A431 cells, HeLa cells

Notable Publications

Author	Pubmed ID	Journal	Application
Lipeng Pei	39901144	Mol Med	WB
Wei Li	39881208	EMBO J	IP
Enyuan Zhang	39012622	J Proteome Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28770-1-AP (KAT2B antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.