

For Research Use Only

Harmonin Polyclonal ANTIBODY



Catalog Number: 51028-2-AP **1 Publications**

Basic Information

Catalog Number: 51028-2-AP	GenBank Accession Number: BC016057	Purification Method: Antigen affinity purification
Size: 150UL	GeneID (NCBI): 10083	
Source: Rabbit	Full Name: Usher syndrome 1C (autosomal recessive, severe)	
Isotype: IgG	Calculated MW: 552 aa, 62 kDa	
Immunogen Catalog Number: AG0513		

Applications

Tested Applications:
ELISA

Cited Applications:
IHC

Species Specificity:
human, mouse

Cited Species:
mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Grimsley-Myers Cynthia M CM	20016102	J Neurosci	IHC

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data